



Approved Mask Cleaning agents

- Steri-7
- Protek EnviroCair
- Equissept
- Virkon
- Dettol
- Accel TB
- Safe4



Approved Med Cup Cleaning agents

- Fairy hand dishwashing up liquid
- Dawn dishwashing liquid
- Dreft dishwashing liquid
- Ecover hand dishwashing up liquid
- Henkel's Pril washing up liquid



We do not recommend any other cleaners

If you have any questions or are unsure if the product you have is suitable for use please contact the **Nortev Customer Service Team** at info@nortev.com or +353 (0) 91 799778 and we will be happy to help.